

# Prior Focus Controller

## Next Generation Motor Controllers for Affordable Focus Control.



The Focus only controller from Prior offers a most affordable solution for microscope focus automation today.

This focus only option is comprehensive in its ability to automate the widest range of microscopes.

### Features:

- Compact Design 160x125x60mm
- Controls:
  - ◇ Manual Focus Knob
  - ◇ Multi Speed Options
  - ◇ Fast Up Buttons
  - ◇ Fast Down Buttons
  - ◇ Position Readout
- Focus Encoder Input
- USB & RS232 Communications
- Compatible with most imaging software packages

Prior motorised stages have a reputation for quality and performance. As an ISO 9001:2000 accredited company Prior equipment is designed and manufactured to the highest quality standards. Prior provides full support and service both direct and indirect through a professional, knowledgeable and extensive global dealer network.

### Options available:

Part Number	Description
ES10ZE	Focus only controller for encoded or un-encoded focus motors with screen, manual focus knob and RS232 and USB communications. Requires PS3H122 style focus motor.

**Prior Scientific Limited**  
Cambridge, UK.  
T. +44 (0) 1223 881711  
E. [uksales@prior.com](mailto:uksales@prior.com)

**Prior Scientific Inc.**  
Rockland, MA, USA.  
T. +1 781-878-8442  
E. [info@prior.com](mailto:info@prior.com)

**Prior Scientific GmbH**  
Jena, Germany.  
T +49 (0)3641 675 650  
E. [jena@prior.com](mailto:jena@prior.com)

VISIT PRIOR ON THE WEB AT [WWW.PRIOR.COM](http://WWW.PRIOR.COM)

